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INDICATIONS FOR USE

The Gemini Line prosthetic instrument set consists of a series of manual instruments, classified as class I medical devices, without measuring purposes, reusable, non-sterile.

DESCRIPTION OF THE DEVICES AND TECHNICAL SPECIFICATIONS

The line includes:

- non-sterile and sterilisable instruments, intended for use in the oral cavity, up to the pharynx, for a timeframe of within 60 minutes and intended for tightening/unscrewing prosthetic screws and positioning/extracting prosthetic components.

The devices are made from materials fit for each of their intended uses, in order to withstand the stress to which they are subjected:

- Invasive instruments:
 - AISI 303: V4 Blocker, MUA Positioner, Analogue Positioner;
 - AISI 420MOD - AISI 303: Dynamometric wrench adapters;
 - AISI 420 MOD: Dynamometric Wrench Drivers, Manual Torque Drivers, MUA Converter Torque Drivers, Sly System Extractor, Sly System Activator.
- Non-invasive instruments:
 - Delrin and AISI 303: V4 Regolo;
 - Resin: Angle Tester.

Further technical details can be found in the New Ancorvis catalogues or on www.newancorvis.eu

CONTRAINDICATIONS

- Do not use the devices for purposes other than their intended use;
- Do not use on patients who are allergic or intolerant to the material constituting the MD.

METHOD OF USE AND WARNINGS

- Before each use, check the integrity of the devices: if they show detached, oxidised, fractured, broken or worn parts, do not use;
- Before use, check the information concerning the compatibility of the devices. The responsibility for checking the compatibility of the devices with prosthetic and/or implant components of other establishments lies with the user (Dentist or Dental Technician, as appropriate);
- Use and handling of devices is reserved for Dental Technicians and Dentists who are properly trained and qualified to practise in premises designated for such activities;
- When tightening the prosthetic components, comply with the torque suggested by the respective manufacturers of the components: Unsuitable torques may compromise the mechanical structure of the components used and/or compromise the stability of the connection;
- The measurement scale available on some of the instruments is purely indicative.

CLEANING, WASHING, AND STERILISATION

Caution: the MDs are supplied non-sterile and must be cleaned, disinfected and sterilised before each use as described below.

The protocol described below is valid only for MDs produced by New Ancorvis. In order to ensure the effectiveness of the process, we therefore recommend processing these MDs without mixing them with MDs from other manufacturers, with other materials, sizes, shapes and intended uses.

- **Functional checks:** before carrying out the procedure described below, check the condition of the MDs: if they show signs of wear and/or damaged parts, dispose of the product.
- **Prewash:** before each use, the MDs must be immersed in a solution of warm water (<40°) and detergent/disinfectant (such as Bactizyme) for 10-15 minutes. Take the MDs and brush them carefully, with a soft bristle nylon brush, under running water in order to remove any organic residues. Follow the detergent manufacturer's instructions to determine the amount and/or dilution of the product. When brushing, pay special attention to the grooves, holes and threads, with the aid, if necessary, of special pipe

cleaners to reach the most difficult parts. Visually check the MDs after the procedure: if there are residues, repeat until they disappear.

- **Ultrasonic cleaning/disinfection:** the MDs cleaned as described, must be placed in an ultrasonic tank or disinfectant and subjected to a cycle of at least 10 minutes at 40° C, using a detergent/enzyme-based disinfectant and/or quaternary ammonium salts (such as Bactizyme) in the dilution suggested by the product manufacturer. Avoid using chlorine-based solutions (e.g. saline solutions). Rinse for at least 1 minute under running water and dry with compressed air or a lint-free cloth.
- **Packaging:** bag the MDs in pouches compatible with the sterilisation cycle described.
- **Sterilisation:** sterilise with saturated steam in an autoclave at 134° for at least 5 minutes.
- **Storage:** the sterilised pouches must be stored in a cool, dry place away from heat sources, light and dust. Open the packages only when using the contents.
- **Number of uses:** the MDs are validated for numerous cycles, depending on their state of wear, at the end of which they must be disposed of (cycle means use, pre-wash, ultrasonic cleaning/disinfection and sterilisation). If the MDs show damage or signs of wear before this limit, they must in any case be disposed of.

DISPOSAL

The devices must be disposed of in compliance with relevant legislation in force: if they are used, they must be disposed of in accordance with biological waste disposal legislation in force.

DISCLAIMER

The Gemini Line prosthetic instruments must be used by specialist personnel who understand clinical use protocols and can recognise any defects or signs of wear on the devices.

New Ancorvis denies all liability for direct and/or indirect damages arising from the user's inexperience, any changes made by the user to the original form of the devices, their improper use and/or incorrect storage and treatment.

ADVISORY CONCERNING SERIOUS INCIDENTS

Should a serious accident occur while using the devices or in the context of their use, the patient and/or the user must report it to the competent Authority in the country where the event occurred and to the Manufacturer NEW ANCORVIS srl at the following address info@newancorvis.eu, specifying the code and batch of the product involved.

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