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Rev. 02\_16.11.2022

## TECHNICAL DATA SHEET

<b>Commercial name</b>	Gemini Line - V4 System - V4 Base																											
<b>Code</b>	Consult the catalogue to view the various configurations available.																											
<b>Intended use</b>	<p>The medical device, which is secured to the superstructure by anaerobic cement, allows a stable and uniform fit (titanium/titanium) between the fixture and prosthesis.</p> <p>The V4 Base, along with its accessories, creates the multi-purpose V4 System</p>																											
<b>Description of the device and how to use it</b>	<p>The MD consists of a tri-lobed cylindrical cannula with a pass-through hole in Titanium Gr.5.</p> <p>The cannula is adjusted in height by the operator using the specific tool.</p> <p>The tri-lobed shaped of the cannula evenly distributes the cement and precisely positions the converter in the prosthesis.</p> <p>The correct management of the cement thickness to use is determined by the use of specific libraries.</p> <p>The pass-through hole holds the screw, used to fix the connector to the underlying implant, which can be either a straight hole method or F.A.S. method (Free Angle System).</p> <p>The ti-base is supplied with a screw available in both standard and High Grip versions.</p>																											
<b>Contact area features and application time</b>	<p>The MD can come into contact with the mucose of the oral cavity and biological fluids.</p> <p>Application time more than 30 days.</p>																											
<b>Classification</b>	IIb (Annex VIII MDR 745/2017)																											
<b>CE Marking</b>	Certificate CE <sub>0425</sub> – No. 0425-MED-003980-00																											
<b>Technical features</b>																												
<b>Material</b>	<p>Titanium Gr. 5 (EN TiAl6V4)</p> <table border="1" style="width: 100%; text-align: center;"> <thead> <tr> <th>C%</th> <th>Fe%</th> <th>O%</th> <th>N%</th> <th>H%</th> <th>Al%</th> <th>V%</th> <th>Y%</th> <th>Ti%</th> </tr> </thead> <tbody> <tr> <td>Max</td> <td>Max</td> <td>Max</td> <td>Max</td> <td>Max</td> <td>5,50</td> <td>3,50</td> <td>Max</td> <td>saldo</td> </tr> <tr> <td>0,08</td> <td>0,25</td> <td>0,13</td> <td>0,05</td> <td>0,012</td> <td>6,50</td> <td>4,50</td> <td>0,005</td> <td></td> </tr> </tbody> </table>	C%	Fe%	O%	N%	H%	Al%	V%	Y%	Ti%	Max	Max	Max	Max	Max	5,50	3,50	Max	saldo	0,08	0,25	0,13	0,05	0,012	6,50	4,50	0,005	
C%	Fe%	O%	N%	H%	Al%	V%	Y%	Ti%																				
Max	Max	Max	Max	Max	5,50	3,50	Max	saldo																				
0,08	0,25	0,13	0,05	0,012	6,50	4,50	0,005																					
<b>Screw tightening</b>	See product label																											
<b>Single use</b>	YES																											
<b>Status</b>	Non-sterile																											
<b>Sterilisation method</b>	In steam autoclave at 134°C for at least 5 minutes, as per validated cycle. The device must be cleaned and washed before sterilisation using a validated method based on the residue of the user processing. Do not use dry or chemical autoclaves.																											
<b>Storage method</b>	Store in a cool, dry place																											
<b>Presence of hazardous substances</b>	The product does not contain SVHC substances. (Substances of Very High Concern). Ref. REACH Regulation - art. 59 - par. 10)																											

Rev. 02\_16.11.2022