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TECHNICAL DATA SHEET

Commercial name	Gemini Line – V4-Ti Base / V4-Ti Base F.A.S. / V4-Ti Base HG/ V4-Ti Base HG F.A.S
Code	V4[x][c] / V4[x][c] / V4HG[x][c] / V4HG[x][c] Legend: [x]= compatibility and dimensional specifications [c]= colour if applicable Consult the catalogue to view the various configurations available.
Intended use	The medical device, which is secured to the superstructure by anaerobic cement, allows a stable and uniform fit (titanium/titanium) between the fixture and prosthesis. The V4 Ti Base, along with its accessories, creates the multi-purpose V4 System. In diatoric restorations, the use of V4+ is recommended.
Description of the device and how to use it	The MD consists of a tri-lobed cylindrical cannula with a pass-through hole in Titanium Gr.23. The cannula is adjusted in height by the operator using the specific tool. The tri-lobed shaped of the cannula evenly distributes the cement and precisely positions the converter in the prosthesis. The correct management of the cement thickness to use is determined by the use of specific libraries. The pass-through hole holds the screw, used to fix the connector to the underlying implant, which can be either a straight hole method or F.A.S. method (Free Angle System). The ti-base is supplied with a screw available in both standard and High Grip versions.
Contact area features and application time	The MD can come into contact with the mucose of the oral cavity and biological fluids. Application time more than 30 days.
Classification	IIb (Annex VIII MDR 745/2017)
CE Marking	Certificate CE 0425 – No. 0425-MDR-030026-02 (exp. 09.05.2028)
Technical features	
Material	Titanium (Ti6Al4V-ELI) Regulatory Standard: ASTM F136
Screw tightening	See product label
Single use	YES
Status	Non-sterile
Sterilisation method	In steam autoclave at 134°C for at least 5 minutes, as per validated cycle. The device must be cleaned and washed before sterilisation using a validated method based on the residue of the user processing. Do not use dry or chemical autoclaves.
Storage method	Store in a cool, dry place

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Presence hazardous substances	of	The product does not contain SVHC substances. (Substances of Very High Concern). Ref. REACH Regulation - art. 59 - par. 10)
--	-----------	---