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Rev. 04_10.05.2023

TECHNICAL DATA SHEET

Commercial name	Gemini Line - Reinforced support V4 LEGO/V4 LEGO BAR.
Code	V4LEGO-[x] Legend: -[x]= compatibility and dimensional specifications Consult the catalogue to view the various configurations available.
Intended use	The Medical Device identifies a series of components which, when combined, based on the prosthetic project to be created, allow a support structure to be assembled for a temporary prosthesis for application on the implants.
Description of the device and how to use it	The components are manufactured using the Laser Melting additive printing technique and are available in either Titanium Gr. 5 or CrCo. It is available in two versions: <ul style="list-style-type: none"> - V4 Lego A: the MD consists of a cylindrical body that can be attached to the V4 Base. Anaerobic cement can be injected through the cylindrical hole to allow it to be attached to the underlying V4 Base. The external supports are intended to accommodate wire/plate/metal bar for the fabrication of the temporary restoration by the dental technician. - V4 Lego Bar: consisting of different types of components, which can be ordered individually or in kits, divided into threaded joints with variable angle - telescopic elements - end joints - cantilever joints.
Contact area and application time	The MD is incorporated by the aesthetic coating of the temporary prosthesis and so, never comes into direct contact with the mucose of the oral cavity. It, therefore, stays in position in the oral cavity for the entire time the temporary prosthesis is in use (from 30 days to 1 year).
Classification	Ila (Annex VIII MDR 745/2017)
CE Marking	Certificate CE ₀₄₂₅ - No. 0425-MED-030026-00 (exp. 09.05.2028)
Technical features	
Material	Titanium Gr. 5 or CoCr
Single use	YES
Status	Unsterilised
Sterilisation method	In steam autoclave at 134°C for at least 5 minutes, as per validated cycle. The device must be cleaned and washed before sterilisation using a validated method based on the residue of the user processing. Do not use dry or chemical autoclaves.
Storage method	Store in a cool, dry place
Presence of hazardous substances	The product does not contain SVHC substances. (Substances of Very High Concern). Ref. REACH Regulation - art. 59 - par. 10)

Rev. 04_10.05.2023