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INDICATIONS FOR USE

The prosthetic components are designed and manufactured to allow a correct and optimal connection of the prosthetic elements with the various implant components already on the market.

V4 BASE - CONVERTERS FOR CEMENT-RETAINED TECHNIQUE

The converters for the cement-retained technique are made of grade 5 titanium both with non-rotational engagement (for single prosthesis) and with rotational engagement (for multiple prostheses, Toronto Implant Bridge); The height of the cementable canula can be modified by the dental technician, as needed, using the dedicated V4 Reg-olo accessory.

The MD is marketed in single packs, accompanied by the respective definitive prosthetic screw (with or without HG treatment), manufactured by New Ancorvis. The types of converters are produced to be used both for temporary and for permanent prostheses.

Technical specifications

- Single use.
- Grade 5 titanium to ensure excellent resistance to the masticatory force.
- Available in various transmucosal heights to better adapt to the aesthetic needs of the patient.
- To connect to the implant, use the through screw included in the package.

HIGH GRIP SCREWS

The High Grip prosthetic screws are characterised by a special surface coating, on the threaded part, which allows reducing the friction due to the systematic micro-vibrations caused by chewing. They are intended only for permanently fixing the prosthesis/implant in the oral cavity.

The prosthetic screws made by New Ancorvis are compatible and therefore can replace the original screws of the implant manufacturer.

Technical specifications:

- Single use.
- The screws are made of grade 5 titanium for excellent resistance to the masticatory force and treated with a coating.
- The tightening torque to be applied to the screw is indicated on the label.

V4 LEGO/V4 LEGO BAR - REINFORCED SUPPORT FOR TEMPORARY PROsthESIS

V4 Lego BAR identifies a series of components which, combined with each other and based on the prosthetic project to be created, allow assembling a support structure for a temporary prosthesis to be applied on implants.

The components are made with the Laser Melting additive printing technique and are available both in grade 5 Titanium and in CrCo and are divided into joints, telescopic bars and cylinders.

Select joints, telescopic bars and cylinders suitable for the type of project to be carried out.

Apply each joint to the respective V4 Base by screwing until the desired height and orientation is achieved.

Fix the joint or cylinder to the V4 Base by injecting the fixing cement through the appropriate accesses or alternatively by laser sealing.

Compose each telescopic bar, inserting the two components one inside the other until the desired length is reached, taking care to maintain adequate proportions for the coupling to give stability to the assembly.

Cut off the excess part of each bar with a separating disc.

Assemble the bars thus obtained on the respective joints and fix the various components using:

- Tightening cement, following the instructions provided by its manufacturer;
- Laser sealing.

The reinforcing structure is ready to be finished with aesthetic coating.

Technical specifications:

- Single use.
- The components are available both in grade 5 titanium and in CoCr alloy to ensure excellent resistance to the masticatory force.
- Components are available in various configurations to accommodate any planned prosthetic reconstruction.

V4 PICK UP

V4 Pick Up is a transfer that can only be used in combination with the V4 Base.

The MDs are identified by screen printing and anodic colouring to facilitate the V4 Pick Up/V4 Base/Implant correspondence.

To ensure correct impression taking, it is recommended to use a magnifying system to verify that the transfer body is screwed all the way into the V4 Base.

For use in the impression coping configuration, the dedicated coping screw must be used for each implant platform; all coping screws have a knurled part to facilitate manual tightening and engagement for a 1.27mm hex driver (New Ancorvis coding: B).

Technical specifications:

- Multi-use.
- V4 Pick Up and V4 Screw Coping made of grade 5 titanium.

V4 SCAN MARKER

The V4 Scan Marker is the Scan Abutment made exclusively for the V4 Base, therefore use is restricted to this type of product only. The coupling between the two components is univocal as the Scan Marker can be inserted in only one position, thanks to the trilobate body of the cement-retained converter, guaranteeing a precise intraoral scan; the insertion phase is facilitated by the two notches (one on the converter, one on the scan marker) which must coincide. It is essential to use the MD exclusively with the new V4 base component, as should it be necessary to adjust its height, it cannot be scanned digitally. Vertical positioning is guaranteed by a click, which certifies the correct positioning in "Z"; if after several uses the "click" is no longer perceptible, we recommend replacing the V4 Scan Marker. It is advisable to keep the part inclined buccally.

Technical specifications:

- Multi-use.
- The V4 San Marker is made in Peek.
- Product to be used only on the V4 Base.
- Use the specific libraries for this product.
- In the case of creation of Implant Bridge, use this product for no more than two implant elements; for larger workings use the AQ scan abutment.

CONTRAINDICATIONS

- Do not use the devices for purposes other than their intended use.
- Do not use on patients with allergies or intolerances to the device's component materials.

METHOD OF USE AND WARNINGS

- Before each use, check the integrity of the devices: if they present anomalies, do not use.
- Before use, check the information concerning the compatibility of the devices.
- The responsibility for checking the compatibility of the New Ancorvis prosthetic components with the prosthetic components and instruments from other manufacturers lies with the user (Orthodontist or Dental Technician as applicable).

- Use and handling of devices is reserved for Dental Technicians and Dentists who are properly trained and qualified to practise in locations designated for such activities.
- The screws, the screw housings and the platforms must NEVER be modified in any way.
- The torque to be applied to the screws is indicated on the relevant label. Unsuitable torques may compromise the mechanical structure of the components used and/or compromise the stability of the connection.
- SINGLE USE DEVICES. The MDs are single-use: any reuse exposes the user to the risk of cross-infection, the appropriate reference symbols are therefore present on the label and the packaging.

CLEANING, WASHING AND STERILISATION FOR SINGLE-USE MDS

The devices are supplied NON-STERILE: correct cleaning, disinfection and sterilisation is essential to reduce risks for the patient and to prevent cross-infections.

The devices are supplied cleaned and washed according to the New Ancorvis protocol, suitably validated.

Equipment and materials needed for the sterilisation process:

Pouches for sterilisation compatible with the method described below - Pouch sealing machine - Steam autoclave.

Packaging and sterilisation

- The devices are supplied NON-STERILE: before using them on the patient, they must be packaged and sterilised.

Packaging:

- Before undergoing the sterilisation cycle, the MD must be removed from its original packaging and packaged in pouches compatible with the sterilisation method described.
- The pouches must be sealed by the user with a validated method and kept in a cool, dry place until sterilised.

Sterilisation:

- Each pouch thus obtained must be placed in a saturated steam autoclave and subjected to a 134°C cycle for at least 5 minutes.
- Storage: the sterilised pouches must be stored in a cool, dry place away from heat sources, light and dust. Open the packages only when using the contents.

CLEANING, WASHING AND STERILISATION FOR MULTI-USE MDS

Caution: the MDs are supplied non-sterile and must be cleaned, disinfected and sterilised before each use as described below.

The protocol described below is valid only for MDs produced by New Ancorvis. In order to ensure the effectiveness of the process, we therefore recommend processing these MDs without mixing them with MDs from other manufacturers, with other materials, sizes, shapes and intended uses.

- Functional checks: before carrying out the procedure described below, check the condition of the MDs: if they show signs of wear and/or damaged parts, dispose of the product.
- Prewash: before each use, the MDs must be immersed in a solution of warm water (<40°) and detergent/disinfectant (such as Bactizyme) for 10-15 minutes. Take the MDs and brush them carefully, with a soft bristle nylon brush, under running water in order to remove any organic residues. Follow the detergent manufacturer's instructions to determine the amount and/or dilution of the product. When brushing, pay special attention to the grooves, holes and threads, with the aid, if necessary, of special pipe cleaners to reach the most difficult parts. Visually check the MDs after the procedure: if there are residues, repeat until they disappear.
- Ultrasonic cleaning/disinfection: the MDs cleaned as described, must be placed in an ultrasonic tank or disinfector and subjected to a cycle of at least 10 minutes at 40°C, using a detergent/enzyme-based disinfectant and/or quaternary ammonium salts (such as Bactizyme) in the dilution suggested by the product manufacturer. Avoid using chlorine-based solutions (e.g. saline solutions). Rinse for at least 1 minute under running water and dry with compressed air or a lint-free cloth.
- Packaging: bag the MDs in pouches compatible with the sterilisation cycle described.
- Sterilisation: sterilise with saturated steam in an autoclave at 134° for at least 5 minutes.
- Storage: the sterilised pouches must be stored in a cool, dry place away from heat sources, light and dust. Open the packages only when using the contents.
- Number of uses: the MDs are validated for numerous cycles, depending on their state of wear, at the end of which they must be disposed of (cycle means use, pre-wash, ultrasonic cleaning/disinfection and sterilisation). If the MDs show damage or signs of wear before this limit, they must in any case be disposed of.
- Do not use dry or chemical autoclaves: their use could seriously compromise the aesthetic and functional structure of the device, making it unsuitable for its intended use.

DISPOSAL

The devices must be disposed of in compliance with relevant legislation in force: if they are used, they must be disposed of in accordance with biological waste disposal legislation in force.

DISCLAIMER

GEMINI RANGE components, made by New Ancorvis, must be used by qualified personnel who are aware of clinical implant and prosthetic restoration protocols and are able to recognise any defects that the devices may present. New Ancorvis denies all liability for direct and/or indirect damages arising from the user's inexperience, any changes made by the user to the original form of the devices, their improper use and/or incorrect storage and treatment.

ADVISORY CONCERNING SERIOUS INCIDENTS

Should a serious accident occur while using the devices or in the context of their use, the patient and/or the user must report it to the competent Authority in the country where the event occurred and to the Manufacturer NEW ANCORVIS srl at the following address info@newancorvis.eu, specifying the code and batch of the product involved.

SSCP (SUMMARY OF SAFETY AND CLINICAL PERFORMANCES)

The document is available on the naloading.it and naspeedyshop.it platforms, which each customer can access with their own exclusive credentials.

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