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Healing Cap for MUA

DEVICE DESCRIPTION AND VARIANTS

Healing Caps are Medical Devices (MDs) designed to protect the MUA (Multi-Unit Abutment) connection over the long term and prevent contamination or alteration prior to the final prosthetic restoration. They are screw-retained, one-piece components available in several configurations to ensure compatibility with different types of MUAs. Please refer to website <https://naspeedyshop.it/index.php?lang=en> to verify the specifications of each product code.

INTENDED USE AND EXPECTED CLINICAL BENEFIT:

These MDs are intended for use during prosthetic treatment in the intermediate phases between surgery, healing, and final restoration placement.

INDICATIONS FOR USE:

The MDs are non-sterile and single-use preformed components intended to be screwed onto the MUA connection to protect it and promote soft tissue healing. They must be removed prior to intermediate try-ins and final prosthetic placement. Application is performed by screwing the device onto the MUA.

CONTRAINDICATIONS:

Do not use in patients allergic or intolerant to the materials composing the devices. Verify the compatibility of the MDs with any existing prosthetic materials in the oral cavity. This MD is not intended to function as an abutment and must not be subjected to masticatory loads.

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**TARGET USERS:**

Use of the MD is restricted to licensed and trained Dental Professionals operating in authorized clinical settings.

TARGET PATIENTS:

The MD is intended for patients with one or more dental implants.

MATERIAL:

CAP: Ti6Al4V-ELI (Titanium alloy)

INSTRUCTIONS FOR USE:

The following procedure is indicative; application must be performed exclusively by qualified healthcare professionals.

In the Dental Laboratory: N/A – Not Applicable

In the Dental Clinic:

1. Sterilize the MD prior to intraoral use according to the specified guidelines (see 'Cleaning, Washing, and Sterilization').
2. After MUA placement, remove any biological residues (blood, tissue) and/or sutures.
3. Screw the MD onto the MUA, ensuring proper seating and, using appropriate magnification tools, verify the marginal seal.

WARNING:

It is the clinician's responsibility to instruct the patient on proper hygiene maintenance.

WORKABILITY LIMITATION:

The MD is non-modifiable.

CLEANING, WASHING AND STERILIZATION

The Medical Devices (MDs) are supplied cleaned and washed according to the validated protocol by New Ancorvis.

The MDs are non-sterile; proper cleaning, disinfection, and sterilization prior to clinical use are essential to minimize patient risk and prevent local infections or inflammation that could compromise the biological integrity of peri-implant tissues.

In the Dental Clinic:

Before application to the patient, any prosthesis containing the MD(s) must be sterilized according to the instructions provided by the dental laboratory.

In cases of intraoral use of a single MD (e.g., during try-ins), follow the procedure below:

- MDs are supplied non-sterile: prior to patient use, they must be packaged and sterilized.
- Before sterilization, remove the MD from its original packaging and place it in sterilization pouches compatible with steam sterilization.
- Pouches must be sealed using a validated method by the user.
- Each pouch must be placed in a saturated steam autoclave and subjected to a 134°C cycle for at least 5 minutes.
- After sterilization, store the pouch in a cool, dry place away from heat, light, and dust. Open the pouch only at the time of use.

COMPATIBILITY

Compatibility with the applied MUA (Brand, Platform Model) is indicated on the MD label.

Responsibility for selecting the appropriate New Ancorvis prosthetic component compatible with implants, prosthetic components, and instruments from other manufacturers lies with the user.

Before use, verify compatibility information, occlusal load, and restoration design to determine the appropriate component.

Prosthetic Instrumentation:

MDs must be used with suitable and compatible prosthetic instruments (e.g., drivers, torque wrenches).

New Ancorvis provides compatible instruments for its MDs: consult the website <https://naspeedyshop.it/index.php?lang=en> for compatibility specifications.

WARNING: All Prosthetic instruments used in the clinic must be sterilized before use, according to the respective Manufacturer's instructions. Do not use worn, damaged instruments or those incompatible with the application procedure.

CAUTIONS:

- The MD must never be modified in any way.
- Improper torque application may compromise the mechanical integrity of the components, connection stability, and sealing.

PRECAUTIONS:

- Do not use MDs outside their intended purpose.
- Before each use, check the integrity of the MDs; if any defects are present, do not use.
- During clinical procedures, pay attention to the small size of the devices to prevent accidental ingestion by the patient.

REUSE AND RELATED PRECAUTIONS:

These devices are single-use only: any reuse exposes the patient to cross-contamination risks.

Therefore, the appropriate symbol is displayed on both the label and the packaging.

STORAGE:

Store in the original packaging until use, in a cool, dry place protected from moisture, dust, or other contaminants.

DISPOSAL:

Dispose of MDs in accordance with applicable regulations. If used, they must be discarded following the regulations for biologically contaminated waste.

WARRANTY:

New Ancorvis provides a lifetime warranty on its compatible prosthetic components.

DISCLAIMER OF LIABILITY:

The Gemini Line components manufactured by New Ancorvis must be used by qualified personnel familiar with implant and prosthetic restoration protocols and capable of identifying any device defects.

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New Ancorvis disclaims any liability for direct and/or indirect damages resulting from user negligence, unauthorized modifications to the original device design, improper use, or incorrect storage and handling.

NOTICE REGARDING SERIOUS INCIDENTS:

If a serious incident occurs during or in relation to the use of the devices, the patient and/or user must report it to the competent authority in the country where the event occurred and to the manufacturer New Ancorvis srl at the following address: qualita@newancorvis.eu, specifying the product code and lot number involved.

SSCP (Summary of Safety and Clinical Performance):

The SSCP document is available on Eudamed (<https://ec.europa.eu/tools/eudamed>). Until the platform is fully operational, it can be requested via email at qualita@newancorvis.eu.

MRI SAFETY INFORMATION (Implantable devices):

The MDs have not been evaluated for safety or compatibility in Magnetic Resonance (MR) environments. They have not been tested for heating, migration, or image artifacts in MR settings. Scanning a patient with an MD may result in injury. **Important:** both the patient and the MR operator must be informed of this risk.

SYMBOLS:

For an explanation of the symbols, please refer to the document 'Legenda simbologia New Ancorvis'.

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