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INDICATIONS FOR USE

The class I Gemini Line prosthetic components are designed and manufactured to enable the creation of a prosthetic design compatible with original implants and components. The Line includes a series of invasive and non-invasive, reusable and non-sterile medical devices.

DESCRIPTION OF THE DEVICES AND TECHNICAL SPECIFICATIONS

Impression transfer: Invasive MDs used in the oral cavity for a timeframe of within 60 minutes and intended to be connected to the implant in order to allow taking the impression with the traditional method. Made in Titanium Gr. 5 or AISI 303 Stainless steel depending on the models;

Transfer screws: Invasive MDs used in the oral cavity for a timeframe of within 60 minutes and intended to fix the transfer to the respective implant. Made in Titanium Gr. 5;

Digital model analogues - Plaster model analogues: Non-invasive MDs intended to replicate the position of the implant on a model obtained by digital detection or on a plaster model. Made in AISI 303 steel or Gr.5 titanium according to the models;

Scan Abutments: the MDs, made of Titanium Gr. 5 can be used:
 - in the oral cavity, when connected to the implant, for a timeframe of within 60 minutes, to detect the position of the implant itself during an intra-oral scan.
 - in the laboratory, when connected with the analogue on the model to detect the position of the implant using optical scanning techniques;

Scan Reply: Non-invasive MDs intended to be fixed by screwing onto the V4 Base connection, positioned inside a trial base, a prototype or an impression in order to detect the three-dimensional position of an implant within the oral cavity, using a Scanner Desktop. Made in Titanium Gr. 5 (if coloured) or in AISI 303 steel;

V4 Scan Marker: devices used in the oral cavity for a timeframe of within 60 minutes and intended to be attached to V4 bonding converters to detect the three-dimensional position of an implant during an intra-oral scan. Made in Peek.

Further technical details can be found in the New Ancorvis catalogues or on www.newancorvis.eu

CONTRAINDICATIONS

- Do not use the devices for purposes other than their intended use;
- Do not use on patients who are allergic or intolerant to the material constituting the MD.

METHOD OF USE AND WARNINGS

- Before each use, check the integrity of the devices: if they show detached, oxidised, fractured, broken or worn parts, do not use;
- Before use, check the information concerning the compatibility of the devices. The responsibility for checking the compatibility of the devices with prosthetic and/or implant components of other establishments lies with the user (Dentist or Dental Technician, as appropriate);
- Use of the V4 Scan Marker is recommended only and exclusively on uncut and/or lowered V4 bonding converters;
- Use and handling of devices is reserved for Dental Technicians and Dentists who are properly trained and qualified to practise in premises designated for such activities.

CLEANING, WASHING, AND STERILISATION

Caution: the MDs are supplied non-sterile and must be cleaned, disinfected and sterilised before each use as described below.

The protocol described below is valid only for MDs produced by New Ancorvis. In order to ensure the effectiveness of the process, we therefore recommend processing these MDs without mixing them with MDs from other manufacturers, with other materials, sizes, shapes and intended uses.

- **Functional checks:** before carrying out the procedure described below, check the condition of the MDs: if they show signs of wear and/or damaged parts, dispose of the product.
- **Prewash:** before each use, the MDs must be immersed in a solution of warm water (<40°) and detergent/disinfectant (such as Bactizyme) for 10-15 minutes. Take the MDs and brush them carefully, with a soft bristle nylon brush, under running water in order to remove any organic residues. Follow the detergent manufacturer's instructions to determine the amount and/or dilution of the product. When brushing, pay special attention to the grooves, holes and threads, with the aid, if necessary, of special pipe cleaners to reach the most difficult parts. Visually check the MDs after the procedure: if there are residues, repeat until they disappear.
- **Ultrasonic cleaning/disinfection:** the MDs cleaned as described, must be placed in an ultrasonic tank or disinfectant and subjected to a cycle of at least 10 minutes at 40° C, using a detergent/enzyme-based disinfectant and/or quaternary ammonium salts (such as Bactizyme) in the dilution suggested by the product manufacturer. Avoid using chlorine-based solutions (e.g. saline solutions). Rinse for at least 1 minute under running water and dry with compressed air or a lint-free cloth.
- **Packaging:** bag the MDs in pouches compatible with the sterilisation cycle described.
- **Sterilisation:** sterilise with saturated steam in an autoclave at 134° for at least 5 minutes.
- **Storage:** the sterilised pouches must be stored in a cool, dry place away from heat sources, light and dust. Open the packages only when using the contents.
- **Number of uses:** the MDs are validated for numerous cycles, depending on their state of wear, at the end of which they must be disposed of (cycle means use, pre-wash, ultrasonic cleaning/disinfection and sterilisation). If the MDs show damage or signs of wear before this limit, they must in any case be disposed of.

DISPOSAL

The devices must be disposed of in compliance with relevant legislation in force: if they are used, they must be disposed of in accordance with biological waste disposal legislation in force.

DISCLAIMER

Gemini line class I prosthetic components created by New Ancorvis must be used by specialist personnel who understand clinical use protocols and can recognise any defects or signs of wear on the devices.

New Ancorvis denies all liability for direct and/or indirect damages arising from the user's inexperience, any changes made by the user to the original form of the devices, their improper use and/or incorrect storage and treatment.

ADVISORY CONCERNING SERIOUS INCIDENTS

Should a serious accident occur while using the devices or in the context of their use, the patient and/or the user must report it to the competent Authority in the country where the event occurred and to the Manufacturer NEW ANCORVIS srl at the following address info@newancorvis.eu, specifying the code and batch of the product involved.

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