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GEMINI LINE – CLASS I DEVICES

GENERAL INFORMATION

All products shown in this document are Medical Devices (MDs) intended for use within a prosthetic restoration project.

The MDs are available in several configurations, depending on the type, including: compatibility with original implants, diameter, scan body, transmucosal portion, threading, and pitch.

Please refer to the website <https://naspeedyshop.it/index.php?lang=en> to check the specifications for each product code.

The MDs are preformed, non-sterile, and reusable.

CONTRAINDICATIONS:

When intended for contact with the patient, do not use in cases of allergies or intolerances to the materials composing the devices.

Verify the compatibility of the medical devices with any existing prosthetic materials already present in the oral cavity.

TARGET USERS

The use of these MDs is restricted to Dental Technicians and Licensed Dentists who are properly trained and authorized to practice in facilities designated for such activities.

TARGET PATIENTS:

These MDs are intended for use in prosthetic projects for patients with one or more dental implants.

WORKABILITY LIMITATION:

The MDs are not modifiable.

COMPATIBILITY: Implants, Prosthetic Components, and Instruments

Before use, verify the compatibility information related to the MDs.

Prosthetic Instruments: The MDs must be used in combination with suitable and compatible prosthetic instruments (e.g., screwdrivers, torque wrenches, etc.).

New Ancorvis provides instruments compatible with its MDs: please refer to the website <https://naspeedyshop.it/index.php?lang=en> to check compatibility specifications.

WARNING:

All prosthetic instruments used in the clinic must be sterilized before use, according to the respective manufacturer's instructions.

Do not use damaged instruments or those incompatible with the application procedure.

CLEANING, WASHING, AND STERILIZATION:

The MDs are supplied clean and washed according to the validated protocol by New Ancorvis. The MDs are non-sterile: proper cleaning, disinfection, and sterilization prior to any intraoral use are essential to reduce patient risk and to prevent local infections and local inflammations that could compromise the biological integrity of peri-implant tissues.

In the Dental Laboratory: During processing phases in the laboratory, the MDs may be cleaned as needed without requiring disinfection and/or sterilization.

New Ancorvis MDs can be disinfected using cold disinfection methods.

In the Dental Clinic: In cases of intraoral use of a single MD (e.g., during try-ins), follow the procedure below:

- The MDs are supplied non-sterile: before use on a patient, they must be packaged and sterilized.
- Before sterilization, the MD must be removed from its original packaging and placed in pouches compatible with steam sterilization method.
- The pouches must be sealed using a validated method by the user.
- Each pouch must then be placed in a steam autoclave and subjected to a 134°C cycle for at least 5 minutes.
- After sterilization, store the pouch in a cool, dry place away from heat sources, light, and dust. Open the pouch only at the time of use.

Reuse and related precautions:

Warning: The MDs are supplied non-sterile and must be cleaned, disinfected, and sterilized before each use on a patient, as described below.

The following protocol is validated exclusively for MDs manufactured by New Ancorvis. To ensure the effectiveness of the process, it is recommended to process only these MDs, without mixing them with devices from other manufacturers, or with different materials, sizes, shapes, or intended uses.

- **Functional Checks:** Before performing the procedure described below, inspect the MDs. If any signs of wear and/or damage are present, the MD must be discarded.
- **Pre-washing:** Before each use, immerse the MDs in a solution of lukewarm water (<40°C) and detergent/disinfectant (e.g., Bactizyme) for 10–15 minutes. Remove the MDs and brush them thoroughly under running water using a soft-bristled nylon brush to eliminate any organic residues. Follow the detergent manufacturer's instructions for dilution and quantity. Pay special attention to grooves, holes, and threads, using appropriate brushes if necessary to reach difficult areas. Visually inspect the MDs after the procedure and repeat if residues are still present.
- **Ultrasonic Cleaning/Disinfection:** Once cleaned, place the MDs in an ultrasonic cleaner or washer-disinfector and run a cycle of at least 10 minutes at 40°C using an enzymatic and/or quaternary ammonium salt-based detergent/disinfectant (e.g., Bactizyme), following the manufacturer's dilution guidelines. Avoid chlorine-based solutions

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(e.g., saline solution). Rinse under running water for at least 1 minute and dry with compressed air or lint-free cloths.

- **Packaging:** Place the MDs in pouches compatible with the described sterilization cycle.
- **Sterilization:** Sterilize in a steam autoclave at 134°C for at least 5 minutes.
- **Storage:** Store the sterilized pouches in a cool, dry place away from heat sources, light, and dust. Open the packaging only at the time of use.
- **Number of Uses:** The MDs are validated for multiple cycles, depending on their wear condition. Each cycle includes use, pre-washing, ultrasonic cleaning/disinfection, and sterilization. If the MDs show damage or wear before reaching the expected number of cycles, they must be discarded.

CAUTIONS:

- MDs must never be modified in any way.
- The torque to be applied to screws (when provided with the MD) is indicated on the MD label. Inadequate tightening torque may compromise the mechanical structure of the components involved, the stability of the connection, and the seal.

PRECAUTIONS:

- Do not use MDs outside their intended purpose.
- Before each use, check the integrity of the MD: if any anomalies, deformations, wear, or damage are present, the MD must be discarded.
- During clinical procedures (if applicable), pay attention to the small size of the devices: prevent accidental ingestion by the patient.

WARNING:

The only thermal treatment allowed for the MD is sterilization.

STORAGE:

Store in the original packaging until use, in a cool, dry place, protected from moisture, dust, or other contaminants.

DISPOSAL:

MDs must be disposed of in accordance with applicable regulations. If used, they must be discarded following the regulations for biologically contaminated waste.

WARRANTY:

New Ancorvis provides a 2-year warranty from the date of purchase.

DISCLAIMER OF LIABILITY:

The Gemini Line components manufactured by New Ancorvis must be used by qualified personnel familiar with implant and prosthetic restoration protocols and capable of identifying any device defects. New Ancorvis disclaims any liability for direct and/or indirect damages resulting from user negligence, unauthorized modifications to the original shape of the devices, improper use, or incorrect storage and handling.

NOTICE REGARDING SERIOUS INCIDENTS:

In the event of a serious incident occurring during or in connection with the use of the devices, the patient and/or user must report it to the competent authority in the country where the event occurred and to the manufacturer NEW ANCORVIS srl at the following address: qualita@newancorvis.eu, specifying the product code and lot number involved.

SYMBOLS:

For an explanation of the symbols, please refer to the document 'Legenda simbologia New Ancorvis'.

DEVICE DESCRIPTION AND VARIANTS – SPECIFIC INSTRUCTIONS FOR USE

SCREWS:

(Basic UDI-DI: 8053876-01-020N5)

These MDs are intended for use within a dental prosthetic restoration project to secure Class I devices.

MATERIAL:

Ti Variant: Ti6Al4V-ELI
AISI Variant: AISI 303

INSTRUCTIONS FOR USE:

Use by screwing the device in place with the aid of dedicated instruments, applying the torque specified on the label.

TRANSFERS:

(Basic UDI-DI: 8053876-01-018NJ)

These MDs are intended to be connected to the implant to enable impression taking using the traditional method.

MATERIAL:

Ti Variant: Ti6Al4V-ELI
AISI Variant: AISI 303

INSTRUCTIONS FOR USE:

Using the screw included in the package, connect the MD to the implant to take the impression.

V4 PICK UP:

(Basic UDI: 8053876-01-018NJ)

These MDs are intended to transfer the spatial position of a dental implant or MUA from the patient's oral cavity (analog workflow) to an analog model.

The MDs function exclusively when applied to the coronal portion of the V4 Ti Base.

Application is performed by screwing the MD onto an intact and unmodified V4 Ti Base before proceeding with any type of impression. The MD is compatible exclusively with the V4 Ti Base.

The responsibility for selecting a V4 Ti Base that is compatible with implants, prosthetic components, and instruments from other manufacturers lies with the user (Dentist or Dental Technician, depending on the case).

Before use, verify all information regarding device compatibility, masticatory load, and restoration design in order to select the appropriate component.

MATERIAL:

Ti6Al4V-ELI

INSTRUCTIONS FOR USE

In the Dental Clinic (Analog Impression):

- 1) Sterilize the MD as described in the section "Cleaning, Washing, and Sterilization".
- 2) The MD must be used only on an intact and unmodified V4 Ti Base. IMPORTANT: The quality of the final restoration is completely compromised if this instruction is not followed.
- 3) Carefully verify the correct connection of the V4 Ti Base to the implant/MUA in the oral cavity before screwing in the V4 PICK UP (radiographic confirmation is highly recommended).
- 4) Screw the MD onto the V4 Ti Base already connected to the implant in the oral cavity, ensuring it is fully seated using the dedicated V4 BLOCKER accessory.
- 5) The MD must not come into contact with adjacent teeth or other V4 PICK UPS to avoid stress, misalignment, tooth displacement, or incomplete scanning of the MD.
- 6) Proceed with the open tray impression technique using the most suitable material, at the clinician's discretion.
- 7) Once the MD is returned to the lab after use, follow the procedures outlined in the "Cleaning, Washing, and Sterilization" section.

WARNING – The Clinician is responsible for:

- Selecting the appropriate clinical procedure for the case;
- Proper cleaning and disinfection of the impression;
- Proper storage and handling of the impression before sending it to the lab.

In the Dental Laboratory (Analog Model Fabrication):

- 1) The Dental Technician is responsible for following the correct procedures for producing the physical model.
- 2) After model fabrication, follow the instructions in the section "Cleaning, Washing – In the Laboratory".

PLASTER ANALOGS:

(Basic UDI-DI: 8053876-01-017NG)

These MDs are replicas of the implant present in the patient. When connected to a Transfer/Coping within an analog impression, they allow the construction of a working model in plaster or resin for the fabrication of implant-supported prosthetic restorations.

MATERIAL:

Ti Variant: Ti6Al4V-ELI

AISI Variant: AISI 303

INSTRUCTIONS FOR USE

In the Dental Laboratory:

- 1) Verify that the MD is clean (free from residues of previous use), and that the connection is intact and clean.
- 2) Ensure that the assembly with the Transfer/Coping is accurate and free of gaps or movement once connected.
- 3) Proceed with the creation of the physical model using the selected material.
- 4) Follow the manufacturer's instructions for the chosen material before proceeding to the next steps.

DIGITAL ANALOGS:

(Basic UDI-DI: 8053876-01-017NG)

These MDs are replicas of the implant present in the patient. When assembled onto a prototyped model (derived from a digital scan file) and using the dedicated libraries, they allow the construction of a working model within a digital workflow for the fabrication of implant-supported prosthetic restorations.

MATERIAL:

AISI Variant: AISI 303

INSTRUCTIONS FOR USE

In the Dental Laboratory:

- 1) Verify that the MD's geometry has not been altered by previous use, and that the connection is intact and clean.
- 2) Assemble the digital analog onto the prototyped model, ensuring correct seating in the printed socket. Use compatible and appropriate tools for engaging the screw.
- 3) Secure the MD from the base of the prototyped model using the supplied screw.

SCANBODY V4 SCAN MARKER:

(Basic UDI-DI: 8053876-01-019NL)

These MDs are intended to transfer the spatial position of a dental implant or MUA from the patient's oral cavity (digital workflow) or from an analog model to CAD design software, enabling CAD-CAM fabrication of screw-retained, cemented, or bonded prosthetic restorations using the corresponding implant libraries.

Application is performed by press-fitting onto the coronal portion of an intact and unmodified V4 Ti Base, proceeding with any type of scanning.

The responsibility for selecting compatible materials with implants, prosthetic components, and other manufacturers lies with the user (Dentist or Dental Technician, depending on the case). Before use, verify the information regarding device compatibility.

MATERIAL:

PEEK

INSTRUCTIONS FOR USE

In the Dental Laboratory (Physical Model Scanning):

- 1) The MD must be used only on an intact, unshortened, and unmodified V4 Ti Base.

IMPORTANT: Failure to follow this instruction will compromise the quality of the final restoration.

- 2) Press-fit the MD onto the V4 Ti Base already screwed into the physical model from the analog or digital workflow.
- 3) The MD material allows scanning without the use of scanning sprays.
- 4) The presence of liquid or solid residues or deformations on the MD will seriously compromise the quality of the process: follow the instructions in the section "Cleaning, Washing – In the Laboratory".
- 5) Perform the scan following the instructions of the desktop scanner used.
- 6) Import the scan into the CAD software and proceed with the design of the restoration(s) using the dedicated library.
- 7) The resulting design can be used for fabrication with the chosen material and production technique.
- 8) Remove the MD from the V4 Ti Base and follow the cleaning instructions in the section "Cleaning, Washing – In the Laboratory".

In the Dental Clinic (Intraoral Scanning):

For restorations on multiple implants, to ensure final scan accuracy, it is strongly recommended to connect all V4 Scan Markers to the corresponding V4 Ti Bases and perform a single scan, strictly following the scanner manufacturer's instructions, except in the case described in point 6b.

- 1) Sterilize the MD as described in the section "Cleaning, Washing, and Sterilization".
- 2) The MD must be used only on an intact, unshortened, and unmodified V4 Ti Base.
IMPORTANT: Failure to follow this instruction will compromise the quality of the final restoration.
- 3) Carefully verify the correct connection of the V4 Ti Base to the implant/MUA in the oral cavity before connecting the V4 Scan Marker (radiographic confirmation is highly recommended).
- 4) Press-fit the MD onto the V4 Ti Base already screwed into the implant in the oral cavity.
- 5) The MD must not contact adjacent teeth or other V4 Scan Markers to avoid stress, misalignment, tooth displacement, or incomplete scanning.
- 6) If occurs:
 - a) Contact with teeth: do not use the MD.
 - b) Contact with other V4 Scan Markers: it is possible to scan each MD individually, strictly following the intraoral scanner's instructions.
- 7) Ensure that the MD does not interfere with the opposing arch during occlusal registration. If interference occurs, remove the MD.
- 8) Once the MD has been removed, proceed with cleaning, washing, and sterilization as described in the section "Cleaning, Washing, and Sterilization".

WARNING: The Clinician is responsible for selecting the appropriate digital procedure for case resolution.

SCANBODY:

(Basic UDI-DI: 8053876-01-019NL)

These MDs are intended to transfer the spatial position of a dental implant or MUA from the patient's oral cavity (digital workflow) (except for Scanbody Mod. SQ) or from an analog model to CAD design software, in order to proceed with the CAD-CAM fabrication of screw-retained, cemented, or bonded prosthetic restorations as they are associated with the corresponding implant libraries.

Application is performed by screwing the MD onto the implant or MUA prior to performing any type of scan.

MATERIAL:

Ti6Al4V-ELI

INSTRUCTIONS FOR USE – In the Dental Laboratory (Physical Model Scanning):

Rev. 12/09/2025

- 1) Screw the MD onto the analog in the physical model from the analog or digital workflow.
- 2) After placement, ensure there are no interferences between the MD and parts of the model.

If interference is present:

- remove excess material from the model causing the interference;
 - repeat the placement.
- 3) Screw the MD onto the analog using the included screw.
 - 4) The MD's surface treatment allows scanning without the use of scanning sprays.
 - 5) The presence of liquid or solid residues or deformations on the MD will seriously compromise the quality of the process.
 - 6) Perform the scan following the desktop scanner manufacturer's instructions.
 - 7) Import the scan into the CAD software and proceed with the design of the restoration(s) using the dedicated library.
 - 8) The resulting design can be used for fabrication with the chosen material and production technique.
 - 9) Remove the MD from the model.

In the Dental Clinic (Intraoral Scanning) (excluding Scanbody Mod. SQ):

For restorations on multiple implants, to ensure final scan accuracy, it is strongly recommended to screw in all Scanbodies and perform a single scan, strictly following the scanner manufacturer's instructions, except in the case described in point 4b.

- 1) Screw the MD onto the implant using the included screw.
- 2) Carefully verify the correct Scanbody/Implant/MUA connection by:
 - visual inspection;
 - radiograph (highly recommended).
- 3) The MD must not contact adjacent teeth or other Scanbodies to avoid stress, misalignment, tooth displacement, or incomplete scanning.
- 4) If occurs:
 - a) Contact with teeth: use a Scanbody of appropriate volume.
 - b) Contact with other Scanbodies: scan each MD individually, strictly following the intraoral scanner's instructions.
- 5) Ensure that the MD does not interfere with the opposing arch during occlusal registration. If interference occurs, remove the MD.
- 6) After removing the MD, proceed with cleaning, washing, and sterilization as described in the section "Cleaning, Washing, and Sterilization".

WARNING: The Clinician is responsible for selecting the appropriate digital procedure for case resolution.

SCANBODY SCAN REPLY
(Basic UDI-DI: 8053876-01-019NL)

Scanbody Scan Reply devices are intended to:

- 1) Be screwed onto the V4 Ti Base connection positioned within a test bar;
- 2) Be screwed onto a transfer present in an impression;
- 3) Be screwed onto a corresponding milled connection.

MATERIAL:
Ti6Al4V-ELI

INSTRUCTIONS FOR USE

In the Dental Laboratory:

If used with a test bar:

- 1) Verify that the MD's geometry has not been altered by previous use and that the connection is intact and clean.
- 2) Check compatibility of the MD with the screw supplied with the prosthetic component or bar.
- 3) Assemble the MD to the test bar, ensuring there are no gaps or movement once connected.
- 4) Use compatible and appropriate tools for engaging the screw.
- 5) Proceed with scanning using a desktop scanner, following the manufacturer's instructions.

If used on a transfer in an impression:

- 1) Verify that the MD's geometry has not been altered by previous use and that the connection is intact and clean.

- 2) Check compatibility of the MD with the Transfer/Coping and the corresponding screw in the impression.
- 3) Assemble the MD to the Transfer/Coping, ensuring there are no gaps or movement once connected.
- 4) Use compatible and appropriate tools for engaging the screw.
- 5) Proceed with scanning using a desktop scanner, following the manufacturer's instructions.

In the Dental Clinic:

If used on a connection in a provisional restoration or pre-existing prosthesis

(**Note:** The success of this procedure depends on the passivity and accuracy of the provisional/pre-existing restoration.)

- 1) Verify that the MD's geometry has not been altered by previous use and that the connection is intact and clean.
- 2) Check compatibility of the MD with the screw supplied with the prosthetic component of the provisional/pre-existing restoration.
- 3) Assemble the MD to the provisional/pre-existing restoration, ensuring there are no gaps or movement once connected.
- 4) Use compatible and appropriate tools for engaging the screw.
- 5) Proceed with extraoral scanning using an intraoral scanner.

Note: The same procedure (up to step 4) can also be performed in the dental laboratory. In this case, use a desktop scanner for scanning.

SCANBODY ONE SHOT:

(Basic UDI-DI: 8053876-01-019NL)

These MDs are intended to transfer the spatial position of a dental implant or MUA from the patient's oral cavity (digital workflow) or from an analog model to CAD design software, in order to proceed with the CAD-CAM fabrication of screw-retained, cemented, or bonded prosthetic restorations as they are associated with the corresponding implant libraries.

They are available in various configurations based on: compatibility with original implants, scan body geometry, and transmucosal height. These MDs are monoblock: the fixation screw is integrated. Due to their geometric characteristics, they are particularly suitable for scanning fully or partially edentulous arches.

MATERIAL:
Ti6Al4V-ELI

INSTRUCTIONS FOR USE

In the Dental Laboratory (Physical Model Scanning):

- 1) Screw the MD onto the analog in the physical model from the analog or digital workflow.
- 2) After placement, ensure there are no interferences between the MD and parts of the model.
 - a. If interference is present:
 - remove excess material from the model causing the interference;
 - repeat the placement.
- 3) Screw the MD onto the analog.
- 4) The MD's surface treatment allows scanning without the use of scanning sprays.
- 5) The presence of liquid or solid residues or deformations on the MD will seriously compromise the quality of the process.
- 6) Perform the scan following the desktop scanner manufacturer's instructions.
- 7) Import the scan into the CAD software and proceed with the design of the restoration(s) using the dedicated library.
- 8) The resulting design can be used for fabrication with the chosen material and production technique.
- 9) Remove the MD from the model.

In the Dental Clinic (Intraoral Scanning):

To ensure final scan accuracy, it is strongly recommended to screw in all Scanbodies and perform a single scan, strictly following the scanner manufacturer's instructions.

- 1) Screw the MD onto the implant.
- 2) Carefully verify the correct Scanbody/Implant/MUA connection by:
 - visual inspection;
 - radiograph (highly recommended).

Rev. 12/09/2025

- 3) The MD must not contact adjacent teeth or other Scanbodies to avoid stress, misalignment, tooth displacement, or incomplete scanning.
- 4) If occurs contact with teeth: use a Scanbody of appropriate volume.
- 5) Ensure that the MD does not interfere with the opposing arch during occlusal registration. If interference occurs, remove the MD.
- 6) After removing the MD, proceed as described in the section "Cleaning, Washing, and Sterilization".

WARNING: The Clinician is responsible for selecting the appropriate digital procedure for case resolution.

SCANBODY SCAN FLAG IQ

(Basic UDI-DI 8053876-01-019NL)

The MDs are intended to digitally transfer the spatial position of multiple dental implants or Multi-Unit Abutments from the patient's oral cavity using an intraoral scanning workflow.

Each individual MD can be used exclusively in combination with the other MDs included in the SCANLOGIQ Kit and its corresponding software.

The MD enables CAD design and subsequent CAD-CAM manufacturing of screw-retained or cemented prosthetic restorations, as it is associated with the corresponding implant libraries.

The devices are available in several configurations, determined by: Compatibility with original implant systems and Shape and geometry of the scan body.

Each MD is supplied with a dedicated fixation screw.

Due to their geometric characteristics, the MDs are suitable for scanning fully edentulous or partially edentulous arches.

The MD is non-sterile: once removed, it must be processed according to the instructions provided in the section "Cleaning, Washing, and Sterilization."

MATERIAL:

Ti6Al4V-ELI

INSTRUCTIONS FOR USE:

For usage procedures, refer to the Instructions for Use of the SCANLOGIQ Kit.

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