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Rev. 03_12/09/2025

TECHNICAL DATA SHEET

Commercial name	Gemini Line – Digital Analoq
Code	Consult the catalogue to view the various configurations available.
Intended use	The Medical Device is intended to replicate the position of the implant, applied to the patient, on 3D-printed resin model. The analogue set in this model is a hollow cylinder and is fastened with a screw using dedicated libraries.
Description of the device and how to use it	The MD consists of a fixed-height cylinder head with a variable diameter; it replicates the implant connection of the corresponding fixture in the upper part while the lower internal thread allows the analogue to be fastened to the prototype of the model using a dedicated screw for this MD. The seat of the prototype is engaging thanks to the lower triangular- (for hexagonal and tri-lobed connection) or square-shaped portion (for octagonal connection).
Contact area features and application time	The MD does not come into contact with the patient.
Classification	I
CE Marking	CE (MDR 2017/745)
Technical features	
Material	AISI 303
Screw tightening	See product label
Single use	No
Status	Non-sterile
Sterilisation method	In steam autoclave at 134°C for at least 5 minutes, as per validated cycle. The device must be cleaned and washed before sterilisation using a validated method based on the residue of the user processing. Do not use dry or chemical autoclaves.
Storage method	Store in a cool, dry place
Presence of hazardous substances	The product does not contain SVHC substances. (Substances of Very High Concern). Ref. REACH Regulation - art. 59 - par. 10)

Rev. 03_12/09/2025