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INTENDED USE, DESCRIPTION OF DEVICES AND TECHNICAL SPECIFICATIONS

The Gemini class I prosthetic component line includes a series of reusable, non-sterile medical devices:

Impression transfer: invasive devices in relation to the oral cavity for a timeframe of within 60 minutes and intended to be connected to the implant in order to allow an impression to be taken using the traditional method. Made from grade 5 titanium and AISI 303 steel;

Transfer screws: invasive devices in relation to the oral cavity for a timeframe of within 60 minutes and intended to secure the transfer to the corresponding implant. Made from grade 5 titanium;

Analogues for digital models - Analogues for plaster models: non-invasive devices designed to replicate the position of a system on a model obtained by digital detection or plaster model. Made from AISI 303 steel or Gr 5 titanium;

Scan Abutment: devices, made from grade 5 titanium can be used in the following cases:

- oral cavity, when connected to the implant, for a timeframe of within 60 minutes, to detect the position of the implant during an intra-oral scan.
- laboratory, when connected with the analogue on the model to detect the position of the system using optical scanning techniques;

V4 Scan Marker: devices used in the oral cavity for a timeframe of within 60 minutes and intended to be attached to V4 bonding converters to detect the three-dimensional position of an implant during an intra-oral scan. Made from Peek;

The devices are made from materials suitable for each of their intended uses to withstand the stresses to which they are subjected.

Further technical details can be found in New Ancorvis catalogues or on www.newancorvis.eu

CONTRAINDICATIONS

- Do not use the devices for purposes other than their intended use;
- Not to be used on patients who are allergic or intolerant to the materials constituting the devices.

WARNINGS AND INSTRUCTIONS FOR USE

- The devices are supplied as **NON STERILE**: invasive devices must be cleaned, disinfected and sterilised before use;
- Use and handling of devices is reserved for Dental Technicians and Dentists who are properly trained and qualified to practise in locations designated for such activities;
- Before each use it is necessary to check the integrity:
 - of the devices: if any parts are detached, oxidised, broken or showing signs of wear, do not use them;
 - of the sterilisation bag (if applicable): if damaged, change the bag and sterilise once more;
- Before use, check the information concerning the compatibility of the devices. The responsibility for checking the compatibility of the devices with prosthetic and/or implant components of other establishments lies with the user (Dentist or Dental Technician, as appropriate);
- Use of the V4 Scan Marker is recommended only and exclusively on uncut and/or lowered V4 bonding converters.

CLEANING, WASHING AND STERILISATION

The devices are supplied **NON STERILE**: proper cleaning, disinfection and sterilisation is essential in order to reduce the risks to the patient and prevent cross-infection. To perform these operations correctly, we recommend following these instructions:

Cleaning and washing

- The devices are supplied clean and washed in accordance with a protocol validated by New Ancorvis;
- After each use, the instruments must undergo a thorough washing and cleaning process;
- Use detergents suited to the selected washing method (mechanical or manual), strictly respecting the dosage recommended by the producer of the detergent;
- After washing, the device must be thoroughly dried and all its parts checked;

Sterilisation

- The devices are supplied **NON-STERILE**: they must be sterilised before each use on the patient. New Ancorvis suggests steam autoclave sterilisation at 134°C for at least 5 minutes (cycle validated by the user);
- After sterilisation, store the devices inside their sealed sterilisation bag until just before their use;
- **Do not use dry or chemical autoclaves:** their use may seriously compromise the aesthetic and functional structure of the device, making it unfit for its intended use.

STORAGE

Store in a cool, dry place away from moisture, dust or other pollutants.

DISPOSAL

The devices must be disposed of in compliance with relevant legislation in force: if they are used, they must be disposed of in accordance with biological waste disposal legislation in force.

DISCLAIMER

Gemini line class I prosthetic components created by New Ancorvis must be used by specialist personnel who understand clinical use protocols and can recognise any defects or signs of wear on the devices.

New Ancorvis denies all liability for direct and/or indirect damages arising from the user's inexperience, any changes made by the user to the original form of the devices, their improper use and/or incorrect storage and treatment.

NOTICE ABOUT SERIOUS ACCIDENTS

If, during use of the devices or in the context of their use, a serious accident should occur, the patient and/or the user must report it to the competent Authority in the country where the event occurred and to the Manufacturer NEW ANCORVIS srl at the following address info@newancorvis.eu, specifying the code and batch of the concerned product.

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