CE 0425

V4 – Ti Base
V4 – Ti Base F.A.S.
V4 – Ti Base HG
V4 – Ti Base HG F.A.S.
V4+ - Ti Base
V4+ - Ti Base F.A.S
V4+ - Ti Base HG
V4+ - Ti Base HG F.A.S

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DEVICE DESCRIPTION AND VARIANTS:

The V4 – Ti Base components are Medical Devices (MDs) designed to connect the prosthetic restoration to the implant interface via bonding/cementation, using the corresponding screw.
They are available in several configurations based on:

- Compatibility with original implant systems
- Abutment height/diameter
- Gingival diameter
- Emergence profile height
- Availability of the F.A.S. screw

These components are intended for use in single-unit implant restorations with anti-rotational (engaging) connections, as well as in multiple-unit implant restorations where rotational (non-engaging) connections are required.
Each MD is supplied with its corresponding prosthetic screw. In the V4 HG Ti Base version, the included screw is coated.

REV. 22/06/2026

Please refer to the website https://naspeedyshop.it/index.php?lang=en to verify the specifications of each-code.

INTENDED USE AND EXPECTED CLINICAL BENEFIT:

These MDs are intended for use in prosthetic dental restorations in partially or fully edentulous patients, to restore masticatory function, phonetics, and aesthetics.

INDICATIONS FOR USE:

These are non-sterile and single-use, preformed components intended to support both temporary and definitive screw-retained prosthetic restorations on implants.
They are connected via screw-retention, enabling integration between the implant and the prosthetic component/structure.

- For single-unit restorations, the engaging version is recommended.
- For multi-unit restorations, the non-engaging version is suggested.
- In diatoric restorations, the use of V4+ is recommended.

CONTRAINDICATIONS:

- Do not use in patients allergic or hypersensitive to the materials composing the devices.
- Verify compatibility with any other prosthetic materials already present in the oral cavity.
- Ensure adequate occlusal/vertical space prior to use.

BASIC UDI: 8053876-01-002N3

TARGET USERS:
Use of these MDs is restricted to Dental Technicians and Dentists who are properly trained and licensed to practice in facilities designated for such activities.

TARGET PATIENT:
These MDs are intended for patients with one or more dental implants.

MATERIALS:

- V4 Ti Base / Screw: Ti6Al4V-ELI
- HG Screw: Ti6Al4V-ELI
- HG Coating: C_2H_2

INSTRUCTIONS FOR USE:
The following procedure is indicative: the use and application of the MD on the patient must be performed exclusively by qualified healthcare professionals.

In the Dental Laboratory:

1. For digital workflows, select the appropriate library for case resolution, either via scan or after converting the scan body; proceed with CAD design.
2. Screw the MD onto the physical model (analog or digital workflow) to verify marginal seal.
3. After placement, ensure there are no interferences between the MD and the model.
 If interferences are present:
 • remove excess material
 • re-seat the MD
4. If necessary, adjust the abutment height using appropriate, well-maintained burs specific to the material type. In digital workflows, follow the library's specified chimney height.
 WARNING: The anatomical portion of the MD must not be modified or treated. Ensure it is smooth, polished, and free of processing residues.
5. Verify marginal seal and the selected fit of the prosthetic restoration on the MD.
6. For intermediate try-ins, temporary cementation is allowed using suitable temporary cements.
7. Proceed with final bonding/ cementation following the cement manufacturer's instructions, considering adhesive and aesthetic outcomes. Remove all excess to avoid clinical failure (e.g., over-contours, marginal leakage).

INNOVA METHOD: This method involves threading the prosthesis to allow screw-retention of the MD, facilitating removal during try-ins or adjustments. Final delivery should still involve bonding/cementation.

WARNING: The only thermal treatment allowed is sterilization.

In the Dental Clinic:

1. Disinfect the complete restoration (integrated with the MD) before intraoral placement, following lab instructions.
2. If intraoral bonding/cementation is performed, sterilize the MD beforehand as per the 'Cleaning, Washing, and Sterilization' section.
3. In both cases, verify proper seating and connection of the MD to the implant.
4. Screw the MD(s) onto the implant(s), adhering to the torque specified on the label to avoid damage.
5. Verify marginal seal and passive fit of the restoration on the implant.
 WARNING: Ensure compatibility of all prosthetic materials used.
6. After final placement, seal the screw access channel.

WARNING:
It is the clinician's responsibility to:

- Select appropriate materials and clinical procedures for sealing the screw channel.
- Instruct the patient on proper hygiene maintenance.

WORKABILITY LIMITATION:
The height of the abutment chimney may be modified by the Dental Technician and/or the Clinician, provided the minimum height defined by

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the available digital libraries is respected, and in any case, not reduced below **3 mm**.

CLEANING, WASHING, AND STERILIZATION:

The MDs are supplied cleaned and washed according to a validated protocol by New Ancorvis.
They are NON-STERILE: proper cleaning, disinfection, and sterilization prior to clinical use are essential to reduce patient risk and to prevent infections and local inflammations that could compromise the biological integrity of peri-implant tissues.

In the Dental Laboratory:
During laboratory processing, MDs may be disinfected using cold disinfection methods.
The prosthesis, manufactured in the laboratory, containing New Ancorvis MDs, must be disinfected before delivery to the dental clinic.
It is the Dental Technician's responsibility to select a disinfection method compatible with all prosthetic materials used.

In the Dental Clinic:
Prior to patient use, the prosthesis containing the MD(s) must be sterilized according to the instructions provided by the laboratory.
In the case of intraoral use of a single MD (e.g., during try-ins), follow this procedure below:

- MDs are supplied NON-STERILE: they must be packaged and sterilized before use on the patient.
- Before sterilization, remove the MD from its original packaging and place it in sterilization pouches compatible with the described steam sterilization method.
- Pouches must be sealed using a validated method by the user.
- Each pouch must be placed in a steam autoclave and subjected to a 134°C cycle for at least 5 minutes.
- After sterilization, store the pouch in a cool, dry place away from heat, light, and dust. Open the pouch only at the time of use.

COMPATIBILITY: Implants, prosthetic components, and instrumentation

The compatibility of the MD with the original implant system (Brand, Model, Platform) is indicated on the MD label.
It is the responsibility of the user (Clinician or Dental Technician) to select the appropriate New Ancorvis prosthetic component compatible with implants, prosthetic components, and instruments from other manufacturers.
Prior to use, verify the information regarding device compatibility, masticatory load, and restoration design to determine the appropriate component.

Prosthetic Instruments:
MDs must be used in combination with suitable and compatible prosthetic instruments (drivers, adapters, torque wrenches).
Additional tools include:

- **V4 Regolo – Class I**: intended for Dental Technicians to reduce the height of the Ti-base chimney according to the selected cad library.
- **V4 Pick Up – Class I**: intended for Clinicians. It is an accessory to be assembled onto the V4 Ti-base (not reduced in height) to convert it into an impression transfer for analogic elastomeric impressions, enabling the accurate transfer of the implant position.
- Ensure that the components are securely connected and that the assembly is stable within the impression material.
- **V4 Blocker – Class I**: intended for Clinicians. Accessory designed to facilitate assembly and provide enhanced control and stabilization of the V4 Pick Up on the Ti-base during impression procedures.
- **V4 Scan Marker – Class I**: intended for both Clinicians and Dental Technicians. It is an accessory to be assembled onto the V4 Ti-base (not reduced in height) to convert it into a scanbody for digital impression acquisition.

 New Ancorvis provides compatible instruments for its MDs. Refer to https://naspeedyshop.it/index.php?lang=en for compatibility specifications.
 All instruments used in the clinic must be sterilized prior to use, according to the respective manufacturer's instructions. Do not use damaged instruments or those incompatible with the application procedure.

CAUTIONS:

- Screws, screw seats, and connections must NEVER be modified.
- Improper torque may compromise the mechanical integrity of the components involved, the stability of the connection, and the seal.

PRECAUTIONS:

- Do not use MDs outside their intended purpose.
- Before choosing the component, carry out an accurate assessment of the dental position where it will be placed in the oral cavity and the masticatory load it will need to withstand.
- Before each use, check the integrity of the MDs: do not use if any anomalies are detected.
- During clinical procedures, pay attention to the small size of the devices: avoid the risk of accidental ingestion by the patient.

REUSE AND RELATED PRECAUTIONS:

These devices are single-use only: any reuse exposes the patient to the risk of cross-contamination.
Therefore, the appropriate symbol is displayed on both the label and the packaging.

STORAGE:

Store in the original packaging until use, in a cool, dry place, protected from moisture, dust, or other contaminants.

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MDs must be disposed in accordance with applicable regulations. If used, they must be discarded following current regulations for biologically contaminated waste.

WARRANTY:

New Ancorvis provides a lifetime warranty on its compatible prosthetic components.

DISCLAIMER OF LIABILITY:

The Gemini Line components by New Ancorvis must be used by qualified professionals familiar with implant and prosthetic protocols and capable of identifying device defects.
New Ancorvis disclaims any liability for direct or indirect damage resulting from user negligence, unauthorized modifications to the original device shape beyond indicated limits, improper use, or incorrect storage and handling.

NOTICE REGARDING SERIOUS INCIDENTS:

If a serious incident occurs during or in relation to the use of the devices, the patient and/or user must report it to the competent authority in the country where the event occurred and to the manufacturer at the following address: qualita@newancorvis.eu, specifying the product code and lot number involved.

SSCP (Summary of Safety and Clinical Performance):

The SSCP document is available on EUDAMED (https://ec.europa.eu/tools/eudamed). Until the platform is fully operational, it can be requested via email at qualita@newancorvis.eu.

MRI SAFETY INFORMATION (For implantable devices):

These MDs have not been evaluated for safety or compatibility in a magnetic resonance (MR) environment.
They have not been tested for heating, migration, or image artifacts.
Scanning a patient with a MD may cause injury.
Important: Both the patient and MR operator must be informed of this risk.

SYMBOLS:

For an explanation of the symbols, please refer to the document 'Legenda simbologia New Ancorvis'

V4 – Ti Base
V4 – Ti Base F.A.S.
V4 – Ti Base HG
V4 – Ti Base HG F.A.S.
V4+ - Ti Base
V4+ - Ti Base F.A.S
V4+ - Ti Base HG
V4+ - Ti Base HG F.A.S

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